


Connection HPT home



Advice for therapists:

- Play the Test-CD via the HPT home- the headphones are plugged in the way shown above.
- turn the microphone knob on 0 (to the left)
 - Select 1. "Highfrequency filtering" and then:
 - turn filtering on 0% (all other values don't have effects then)
 - conduct the volume turning as described in the Test-CD
 - mark the position of the switch for volume CD and headphones (for future repetition)


Adjustment of HPT home (Training with Level 1-6)

Switching your on-off switch to ON.

1. Wait, till you see 'Please select' on the display:


2. choose '3. Sequencing'

With this option no more adjustments for the programme have to be made. If you want to use the device more individually, choose the 2. Lateral Hf.-Filt. option. In this mode you have to adjust the values step by step. This is only useful, when you have been provided with the according tables for the values by your supervisor/therapist.



3. Chose level one for week 1 and then keep changing to level 6 for week (upper +\- button)

3.1. **When you are able to conduct the training 3 times a day,** you can choose level one for morning, afternoon and evening .. up to week 6 morning, afternoon, evening.




3.2. **When you can only conduct the training 2 times a day,** choose mornings level 1 and for the afternoon session choose level 1 full. Up to week 6 where again you choose level 6 full during the afternoon session.

3.3. **When you can only do the training 1 a day,** choose level 1 full for the first week.. Till level 6 complete for the 6th week.

4. Start the training by pressing the + key on the lower right hand side

4.1. a little arrow which is slowly moving to the right will appear under the picture

4.2. Pause the training by pressing the menu key to the right. When you press this key again, the training will continue.



5. When the time is over the training will automatically end.

5.1. You can also manually end the training by pressing the menu key to the left (< below the display).


5.2. When pressing the key to the left on the menu again (< below the display) you will end up in the beginning with all the options again.

Additional Information:

Level 1 to 6 in full mode last 30 min, whereas level 1 to 6 for morning, afternoon and evening last 10 minutes each (and are parts of the complete level).

When you want to end the training before the official end, this is OK. It seems to be better to stop during the training then force the participant to keep on going. Try to repeat the same training at a later point in time. We recommend you to be training for at least 30 min a day.

During weeks 1-6 the training consists of listening to music (Baroque- and Mozart music). From week 7 on the speech training will start if indicated. The training with listening to music should still be done 1x a day for 10-20 min (keep level 6 for all the following training weeks). The chosen speech material is used with the microphone. For the speech training one starts again with level 1 (during the 1st week) till level 4 (during week 4) and remain with on level 4 (for week 5 and 6). For speech you should not train in any higher levels than level 4.



Frequently asked questions:

1. Switching not possible:

When the HPT will start up with a given level from 1 to 6 and cannot be changed, the device is in the so-called "locked mode". To get out of this mode:


Keep the lowest "+" key pressed, while you are restarting the device.
You will see the text "The locked mode is deactivated" in the display.

2. Signal fails shortly:

When the volume suddenly decreases and there is shortly no sound while you are talking into the microphone or when listening to CD-parts which are noisy, this is the built-in noise-protector. The noise-protector protects the person listening from too loud noise in the headphone.

To prevent this, turn the switch "microphone" (lower left side) only up to highest 6 and increase the volume headphone up to the necessary value to enable you to hear the microphone voice clearly without any failures.

When the noise-protector is active you will be able to see a "!" sign in the display. You can check it out this way.



3. Original factory settings:

When point 1 and 2 do not show any success, you can **factory reset** the HPT.

Press the second "-" button from the top, while restarting the device, until you see the text "HPT initialised!" in the display.

4. Fixing the device:

Before sending a HPT for repair, you should email/call AUDIVA and describe the faults. Many times we were able to help customers with some easy instructions to get rid of faults.